

IMPORTANT: Please read carefully and keep this information for future use.

This fitting guide is intended for the eye care practitioner, but should be made available to the patient upon request. The eye care practitioner should provide the patient with the wearer's guide that pertains to the patients prescribed lens.

DAILY DISPOSABLE (Hioxifilcon A) COLOR SOFT CONTACT LENS

clearcolor™ 1-day
and
Aryan™ color 1-day

CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED PRACTITIONER.

clearlab™

Trademarks used in this document are owned by Clearlab SG Pte. Ltd. Except otherwise stated, all the trademarks indicated as registered in this document are registered in Singapore and may also be registered in other jurisdictions.

Table of Contents

Description of Lens
Material characteristics
Actions
Indications
Contraindications
Caution
Warnings
Precautions
Adverse reactions
Patient selection
Fitting procedure
Clinical Assessment
Follow-up examinations
Lens Handling (in-office cleaning)
Recommended wearing schedule
Replacement Schedule
Emergencies
Reporting of adverse reactions
How supplied

DESCRIPTION OF LENS

This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is available as a single vision spherical lens. The hydrophilic nature of this material allows the lens to become soft and pliable when immersed in an aqueous solution. This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is fabricated from a nonionic polymer.

MATERIAL CHARACTERISTICS

This Soft Hydrophilic Color Soft Contact Lens is made of a nonionic lens material, Hioxifilcon A, which is a random co-polymer of 2-Hydroxyethyl methacrylate (2-HEMA) and glycerol methacrylate (HEMA/GMA) cross-linked with ethylene glycol dimethacrylate. The co-polymer consists of 42% Hioxifilcon A and 58% water by weight when immersed in normal buffered saline solution. The Hioxifilcon A name has been adopted by the United States Adopted Names Council (USAN).

The lenses are available tinted. Lenses are tinted with one or a combination of one or more of the following pigments, 'listed' color additives: C.I. Reactive black 5, C.I. Vat orange 5, Iron oxides, C.I. Pigment green 7, C.I. Vat brown 1, C.I. Vat yellow 3, C.I. Vat blue 6, C.I. vat orange 1, C.I. Vat green 1, C.I. Pigment blue 36, C.I. Pigment violet 23, D&C Green No.6, phthalocyanato (2) copper, D&C Yellow No.10, D&C Red No.17 and Titanium dioxide. Lenses that contain a unique tinting pattern are subsequently processed to incorporate the 'listed' color additives, and contain only the amount of color additive needed to accomplish the intended cosmetic effect.

In the hydrated state, the lens conforms to the curvature of the eye covering the cornea and extending slightly beyond the limbus forming a transparent or tinted optical surface. The hydrophilic properties of the lens require that it be maintained in a fully hydrated state in a solution compatible with the eye. If the lens dries out, it will become hard and appear somewhat warped however, it will return to its proper configuration when completely rehydrated in the proper storage solution.

Chord Diameter:	12.8 mm to 15.0 mm
Center Thickness:	.080 mm to .266 mm
Base Curve:	8.0 mm to 9.4 mm
Spherical Powers (spherical lens)	-10.00 Diopters to +6.00 Diopters

The hydrophilic characteristics allow aqueous solutions to enter the lens and in its fully hydrated state the lens is approximately 58% water by weight. The physical properties of the lens are:

Refractive Index	1.4011 (hydrated)
Light Transmission (tinted)	greater than 95%
Water Content	58%±2%
Oxygen Permeability	25.38 X 10 ⁻¹¹ (cm ² /sec) (ml O ₂ /ml x mm Hg @ 35°C), (revised Fatt method).

ACTIONS

In its hydrated state, this **Soft Hydrophilic Daily Disposable Color Soft Contact Lens**, when placed on the cornea, act as a refracting medium to focus light rays on the retina.

INDICATIONS

This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** (except plano lens) is indicated for daily wear single use only for the optical correction of refractive ametropia (myopia and hyperopia) in phakic or aphakic persons with non-diseased eyes who may have 1.00D or less of astigmatism.

The lens is intended to be worn once and then discarded at the end of each wearing period on a daily basis. The patient should be instructed to start the next wearing period with a new lens.

CONTRAINDICATIONS (REASONS NOT TO USE)

DO NOT USE this **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** when any of the following conditions are present:

- Acute and subacute inflammation or infection of the anterior chamber of the eye.
- Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids.
- Severe insufficiency of lacrimal secretion (dry eyes).
- Corneal hypoesthesia (reduced corneal sensitivity), if not-aphakic.
- Any systemic disease that may affect the eye or be exaggerated by wearing contact lens.
- Allergic reactions of ocular surfaces or adnexa that may be induced or exaggerated by wearing contact lens or use of contact lens solutions.
- Any active corneal infection (bacterial, fungi, or viral)
- If eyes become red or irritated.
- Patients unable to follow the daily disposable lens care schedule.
- Advise patient not to wear this **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** while sleeping.

CAUTION

Due to the small number of patients enrolled in clinical investigation of lens, all refractive powers, design configurations, or lens parameters available in the lens material were not evaluated in significant numbers. Consequently, when selecting an appropriate lens design and parameters, the eye care practitioner should consider all characteristics of the lens that can affect lens performance and ocular health, including oxygen permeability, wettability, central and peripheral thickness, and optic zone diameter.

The potential impact of these factors on the patient's ocular health must be carefully weighed against the patient's need for refractive correction. Therefore, the continuing ocular health of the patient and lens performance on the eye should be carefully monitored by the prescribing eye care practitioner.

WARNINGS

Please reference Warning in the Package Insert included at the end of this Fitting Guide.

PRECAUTIONS

Please reference Precautions in the Package Insert included at the end of this Fitting Guide.

ADVERSE REACTIONS

Please reference Adverse Reactions in the Package Insert Included at the end of this Fitting Guide.

PATIENT SELECTION

Patient communication is vital. Patients who require visual correction but cannot adhere to the recommended wearing schedule of this **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** should not be provided with this lens. All necessary precautions and warnings should be discussed and understood by the patient (*Review Package Insert with patient*).

Fitting procedure for this Soft Hydrophilic Daily Disposable Color Soft Contact Lens

FITTING PROCEDURE

1. Pre-fitting Examination
2. Initial lens power selection
3. Initial lens diameter and base curve selection
4. Initial lens evaluation
5. Follow-up care

1. Pre-fitting Examination

A pre-fitting patient history and examination are necessary to:

- Determine whether a patient is a suitable candidate for daily wear contact lens (refer to contraindications)
- Collect and record baseline clinical information to which post-fitting examination results can be compared
- Make ocular measurements for initial contact lens parameter selection

2. Initial Lens Power Selection

- a) Convert the spectacle Rx to minus cylinder forms
- b) Compensate the spectacle Rx for vertex distance if the power is greater than + or - 4.00 diopters
- c) Drop the cylinder
- d) Add + 0.25 diopters to compensate for minus tear lens
- e) If refractive astigmatism exceeds 0.75 diopters, determine equivalent sphere and then compensate for power by adding +0.25 diopters for minus tear lens

3. Initial Lens Diameter and Base Curve Selection

The lens is currently offered in one diameter (14.20mm) and one base curve (8.6)

4. Initial Lens Evaluation

- a) Check Lens Centration, Movement, and Size

The criteria for a well fit lens is one which centers easily after a blink, bridges the limbus and extends onto the sclera about 1.5 millimeters, lags downward about 1 to 2 millimeters on upward gaze, and does not move excessively as a result of blinking or exaggerated eye movements.

After the trial lens settled on the eye (5 - 10 minutes), manipulate the lens using lid pressure and observe for indications of excessive tightness. The lens should move freely and easily with the slightest pressure and return to the centered position when released.

Movement of the lens on the eye is very important in assessing the fit and performance of the lens. In primary gaze, slight vertical post-blinking lens movement should occur. On upward gaze, the lens should sag approximately 1 - 2 millimeters.

- b) Refract Over the Lens and Determine Visual Acuity

Allow approximately 10 minutes for fluid equilibration and patient adaptation prior to over refracting. Determine best visual acuity when final over refraction has been achieved. If good visual acuity cannot be obtained through the lens with spherocylindrical over refraction, re-evaluation of the physical fit should be considered. Trial lens procedure should be repeated with lenses of different base curves.

- c) Determine the Optical Power for the Lens Selected

When the proper physical fit has been determined, convert the over refraction through the diagnostic lens to equivalent sphere and add this to the power of the trial lens. This will provide the final power of the lens.

5. Follow-up Care

- a) Follow-up examinations, as recommended by the eye care practitioner, are necessary to ensure continued successful contact lens wear.
- b) Prior to a follow-up examination, the contact lenses should be worn for at least one continuous hour and the patient should be asked to identify any problems which might be occurring related to contact lens wear.
- c) With lenses in place on the eyes, evaluate fitting performance to assure that **CRITERIA OF A WELL FITTED LENS** continues to be satisfied. Examine the lenses closely for surface deposition and/or damage.
- d) After the lens removal, conduct a thorough biomicroscopy examination.
 1. The presence of vertical corneal striae in the posterior central cornea and/or cornea neovascularization is indicative of excessive corneal edema.
 2. The presence of corneal staining and/or limbal-conjunctival hyperemia can be indicative of an unclean lens, a reaction to solution preservatives, excessive lens wear, and/or a poorly fitting lens.
 3. Papillary conjunctival changes may be indicative of an unclean and/or damaged lens.

If any of the above observations are judged abnormal, various professional judgments are necessary to alleviate the problem and restore the eye to optimal conditions. If the **CRITERIA OF A WELL FITTED LENS** is not satisfied during any follow-up examinations, the patient should be re-fitted with a more appropriate lens.

FOLLOW-UP EXAMINATIONS:

- Within one week of lens dispensing
- After three weeks of lens wear
- After seven weeks of lens wear
- After each six month period of lens wear.

At the follow-up examinations, the patient should report good subjective quality of vision. Adaptation to vision with this **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** should occur almost immediately and should definitely be reported within the first (1 week) follow-up visit. At these follow-up visits the practitioner should:

1. Check distance and near acuity with lens in place.
2. Over-refract to verify lens prescription.
3. Observe the position of the lens on the cornea. The lens should be centered and move on upward gaze and with a blink.
4. Evert the lids to examine the tarsal conjunctiva and check for incidence of giant papillary conjunctivitis.
5. Remove the lens. Check corneal curvature. There should be no substantial changes in either meridian.

6. Perform a slit-lamp examination with and without Fluorescein. Check for corneal edema, corneal abrasion, vascularization, corneal infiltrates, and perilimbal injection. Reinsert the lens only after all residual Fluorescein has dissipated from the eye.
7. Clean the lens with a prophylactic surfactant cleaner, and examine for deposits, foreign bodies or physical imperfections of the lens surface.

LENS HANDLING (in-office cleaning)

Wash and rinse hands thoroughly, making certain all soap residues have been rinsed away before drying with a lint-free towel. It is suggested to wet the lens while in the eye using wetting drops before removal. Always start with the right lens first in order to avoid mixing the lens. In removing the lens, try to avoid touching the inside (concave) surface of the lens. It is possible, though not likely, that the lens might be inside out; therefore, check the lens by placing it on the index finger and examine its profile. If the edges of the lens tend to point outward, the lens is inside out. After removing the lens from its container assure that it is clean, clear and wet.

Each **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is received in the eye care practitioner's office in a sterile blister pack. To assure sterility, the blister pack should not be opened until ready for use.

To open the blister pack, pull back the lid where indicated. Upon removing the top cover of the blister pack, the lens may be removed and is ready for use.

This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is not reused in diagnostic procedures.

CLEANING

This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is designed as a Daily Disposable lens.

The lens is intended to be worn once and then discarded at the end of each wearing period. The patient should be instructed to start the next wearing period with a new lens. Emergency lens cleaning and disinfection is not recommended. The patient should be reminded to have replacement lenses or back-up spectacles available at all times.

RECOMMENDED WEARING SCHEDULE

Close professional supervision is recommended to ensure safe and successful contact lens wear. If the patient complains of discomfort, decreased vision, ocular injection or corneal edema, the lens should be removed and the patient scheduled for examination. The problem may be relieved by putting the patient on a different wearing schedule or possibly by refitting the lens.

Patients tend to overwear the lens initially. It is important not to exceed the initial wearing schedule. Regular check-ups, as determined by the eye care practitioner, are also extremely important. The maximum suggested wearing schedule for the

This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is suggested below.

<u>DAY</u>	<u>HOURS</u>
1	6
2	8
3	10
4	12
5	14
6	All Waking hours*

STUDIES HAVE NOT BEEN COMPLETED TO SHOW THAT THIS Soft Hydrophilic Daily Disposable Color Soft Contact Lens IS SAFE TO WEAR DURING SLEEP.

This **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** is intended to be worn once and then discarded at the end of each wearing period. The patient should be instructed to start the next wearing period with a new lens.

EMERGENCIES

Emergency lens cleaning and disinfection is not recommended. The patient should be reminded to have replacement lenses or back-up spectacles available at all times.

The patient should be informed that if chemicals of any kind (household products, gardening solutions, laboratory chemicals, etc.) are splashed into the eyes, the patient should:

FLUSH EYES IMMEDIATELY WITH TAP WATER AND IMMEDIATELY CONTACT THE EYECARE PRACTITIONER OR VISIT A HOSPITAL EMERGENCY ROOM WITHOUT DELAY.

REPORTING OF ADVERSE REACTIONS

Practitioners should report any adverse reactions to this **Soft Hydrophilic Daily Disposable Color Soft Contact Lens** within 5 days to Clearlab SG Pte. Ltd. Additional Package Inserts and Patient Instruction/ Wearer's Guides are available from:

Clearlab SG Pte. Ltd.
139 Joo Seng Road,
Singapore 368362
Tel: +65 6749 1090
Fax: +65 6282 3953
Email: Regulatory@clearlab.com
Website: www.clearlab.com

HOW SUPPLIED

Each lens is supplied sterile in sealed blister packs containing buffered normal saline solution. The blister pack is labeled with the base curve, diameter, dioptric power, manufacturing lot number, and expiration date of the lens. The blister pack is also marked as 'NOT FOR INDIVIDUAL RESALE'.

CAUTION: Federal law restricts this device to sale by or on the order of a licensed practitioner.

clearlab™

Clearlab SG Pte. Ltd.

Printed in Singapore

Revision Date: Dec 2015

Doc Number: S-ASP-070-A

Version Number: 04

© 2015 CLEARLAB SG PTE. LTD. All Rights Reserved.